

ENTHERMICS Medical Systems

Blanket Warming Cabinet



DC150

DC250

DC350

DC400

DC750

120V

OPERATION AND CARE MANUAL



Corporate Headquarters:

W164 N9221 Water Street • P.O. Box 443

Menomonee Falls, Wisconsin 53052-0443 • U.S.A.

PHONE: 262.251.8356
www.enthermics.com

800.T.O.B.WARM U.S.A./CANADA
800.862.9276 U.S.A./CANADA

FAX: 262.251.7067
800.329.8744

TABLE OF CONTENTS

Environmental Conditions	1	Troubleshooting Guide	13
Delivery	1	Service	
Unpacking	1	Parts Lists and Drawings	14-17
Safety Procedures and Precautions	2	Zone Heating Pad Locations	18
Installation		Fuse Replacement	18
General Information	3	Options & Accessories	18
Dimension Drawings	4	Wire Diagram	
Preparation	5	Always refer to wire diagram(s) included with the unit for most current version.	
Electrical Information, Weights & Capacities	5	Warranty	
General Warnings	6	Transportation Damage and Claims	Back Cover
Guidance and Manufacturer's Declaration	6-7	Limited Warranty	Back Cover
Operating Instructions			
Blanket Control Features	8		
Optional Blanket Control Features	9		
Operation Procedures	10		
Care and Cleaning			
Cleaning and Preventative Maintenance			
Protecting Stainless Steel Surfaces	11		
Cleaning Agents	11		
Cleaning Materials	11		
Annual Preventative Maintenance	11		
Clean the Unit Regularly	12		

ENVIRONMENTAL CONDITIONS

Transport and Storage Environmental Conditions (not to exceed 15 days)

- Ambient temperature range of -40° to +70°C (-40° to +159°F)
- Relative humidity range of 10% to 100%, including condensation
- Atmospheric pressure range of 50KPa to 106KPa

Operational Environmental Conditions

- Unit must be allowed to acclimate for a period of at least 24 hours in the environment it will be placed prior to use
- Recommended environmental temperature range is 15°C to 32°C (60°F to 90°F)
- Recommended relative humidity is above 20%, non-condensation

UNPACKING AND SET-UP

DELIVERY

The warming cabinet has been thoroughly tested and inspected to insure only the highest quality unit is provided. Upon receipt, check for any possible shipping damage and report it at once to the delivering carrier. See *Transportation Damage and Claims* section located in this manual.

This appliance, complete with unattached items and accessories, may have been delivered in one or more packages. Check to ensure that all standard items and options have been received with each model as ordered.



Save all the information and instructions packed with the appliance. Complete and return the warranty card to the factory as soon as possible to assure prompt service in the event of a warranty parts and labor claim.

This manual must be read and understood by all people using or installing the equipment model. Contact the service department if you have any questions concerning installation, operation, or maintenance.

NOTE: All claims for warranty must include the full model number and serial number of the unit.

UNPACKING

1. Carefully remove the appliance from the carton or crate.



NOTE: Do not discard the carton and other packaging material until you have inspected the unit for hidden damage and tested it for proper operation.

2. Read all instructions in this manual carefully before initiating the installation of this appliance.

DO NOT DISCARD THIS MANUAL.

This manual is considered to be part of the appliance and is to be provided to the owner or manager of the business or to the person responsible for training operators. *Additional manuals are available from the service department.*

3. Remove all protective plastic film, packaging materials, and accessories from the appliance before connecting electrical power.

CAUTION



TO PREVENT PERSONAL INJURY,
USE CAUTION WHEN MOVING OR
LEVELING THIS APPLIANCE.



WARNING

TRANSPORT SHALL ONLY BE DONE
WITH THE DOORS CLOSED

SAFETY PROCEDURES AND PRECAUTIONS

Knowledge of proper procedures is essential to the safe operation of electrically energized equipment. In accordance with generally accepted product safety labeling guidelines for potential hazards, the following signal words and symbols may be used throughout this manual.



DANGER

Used to indicate the presence of a hazard that **will** cause severe personal injury, death, or substantial property damage if the warning included with this symbol is ignored.



WARNING

Used to indicate the presence of a hazard that **can** cause personal injury, possible death, or major property damage if the warning included with this symbol is ignored.



CAUTION

Used to indicate the presence of a hazard that can or will cause minor or moderate personal injury or property damage if the warning included with this symbol is ignored.

CAUTION

Used to indicate the presence of a hazard that can or will cause minor personal injury, property damage, or a potential unsafe practice if the warning included with this symbol is ignored.

NOTE:


Used to notify personnel of installation, operation, or maintenance information that is important but not hazard related.

1. This blanket warming cabinet is intended for warming cotton blankets **ONLY**. No other use for this device is authorized or recommended.
2. This device is intended for use in commercial establishments where all operators are familiar with the purpose, limitations, and associated hazards of this device. Operating instructions and warnings must be read and understood by all operators and users.
3. Any troubleshooting guides, component views, and parts lists included in this manual are for general reference only and are intended for use by qualified technical personnel.
4. This manual should be considered a permanent part of this device. This manual and all supplied instructions, diagrams, schematics, parts lists, notices, and labels must remain with the device if the item is sold or moved to another location.

NOTE

This unit should not be left unattended for periods of more than 24 hours. In case of absences longer than 24 hours, disconnect the warmer from its power source.

NOTE



For equipment delivered for use in any location regulated by the following directive:

DO NOT DISPOSE OF ELECTRICAL OR ELECTRONIC EQUIPMENT WITH OTHER MUNICIPAL WASTE.

GENERAL INFORMATION

Specifications:

- Single-chamber warming cabinet
- White epoxy-coated steel exterior casing and interior insert
- Single pane, energy efficient, e-coated glass window in door allows for inventory observation
- Easy, hands-free, push-button door release
- Door is fully gasketed and hinged on the right side of the unit
- WarmRight® incorporates a multiple zone warming technology that heats where and when it is needed. All chamber surface temperatures are monitored, providing an efficient balance of heat, low energy consumption and minimal heat loss.
- DC350 & DC750 include a heated center shelf
- Furnished with four (4) 3/4" (19mm) non-skid rubber feet



Control:

- Adjustable temperature range of 37 - 82°C (98 - 180°F)
- Operates in Celsius or Fahrenheit
- Four digit L.E.D. display
- On/Off key
- Up and down arrow keys
- Actual temperature key
- Interior light key
- Built-in speaker for audible feedback
- Built-in control button lock-out feature

Additional features:

- L.E.D. interior lighting casts a comforting blue glow with two (2) different intensity settings, and off mode.
- In the event of a power failure the cabinet will remember its programming and begin to operate as before when power is restored.
- A warming shut-off system, separate from the electronic control, prevents overheating.
- Access point and plug has been added in the back panel to allow the appliance owner to add their data logger or temperature management hardware.

DIMENSIONS



Optional 3" (76mm) casters increases height by 3.6" (91.4mm)

PREPARATION

Before operating the cabinet, clean both the interior and exterior of the unit with a damp cloth and mild soap solution. Wipe with an appropriate disinfectant. Wipe dry with a clean cloth or air dry.

ELECTRICAL INFORMATION & CAPACITIES



The power specifications are located on the unit identification rating tag. This tag is permanently attached to the unit and must be located to verify power requirements.

UL FILE # <input style="width: 80%;" type="text"/>	MADE IN USA
MODEL <input style="width: 80%;" type="text"/>	
WATTS <input style="width: 80%;" type="text"/>	
VOLTS <input style="width: 20%;" type="text"/> 1 PH <input style="width: 20%;" type="text"/> Hz	
SERIAL# <input style="width: 80%;" type="text"/>	

DC150 POWER REQUIREMENTS

120 V.A.C. — 60 Hz, 1 ph
0.6 kW, 5.0 Amps
Safety Class I Equipment
No Applied Parts
Mode of Operation: Continuous



NEMA 5-15P
15A - 125V Plug
Hospital Grade



IPX-0

DC250 POWER REQUIREMENTS

120 V.A.C. — 60 Hz, 1 ph
0.6 kW, 5.0 Amps
Safety Class I Equipment
No Applied Parts
Mode of Operation: Continuous



NEMA 5-15P
15A - 125V Plug
Hospital Grade



IPX-0

DC350 POWER REQUIREMENTS

120 V.A.C. — 60 Hz, 1 ph
0.8 kW, 6.7 Amps
Safety Class I Equipment
No Applied Parts
Mode of Operation: Continuous



NEMA 5-15P
15A - 125V Plug
Hospital Grade



IPX-0

DC400 POWER REQUIREMENTS

120 V.A.C. — 60 Hz, 1 ph
0.8 kW, 6.7 Amps
Safety Class I Equipment
No Applied Parts
Mode of Operation: Continuous



NEMA 5-15P
15A - 125V Plug
Hospital Grade



IPX-0

DC750 POWER REQUIREMENTS

120 V.A.C. — 60 Hz, 1 ph
1.4 kW, 11.7 Amps
Safety Class I Equipment
No Applied Parts
Mode of Operation: Continuous



NEMA 5-20P
20A - 125V Plug
Hospital Grade



IPX-0

Wire diagram is located under top cover of unit.

Clearance requirements:

2" from rear
1" from sides
3/4" from bottom

Grounding reliability can only be achieved when equipment is connected to an equivalent receptacle marked "Hospital Grade."



Protective
Earth
Ground Symbol

Medical Equipment classified by Underwriters Laboratories with Respect to Electrical Shock, Fire and Mechanical Hazards only, in Accordance with UL 60601-1 and CAN/CSA C22.2 No. 601.1.



UL File No.
E201645



Hazardous
Voltage
Present



DANGER



**ENSURE POWER SOURCE
MATCHES VOLTAGE IDENTIFIED
ON APPLIANCE RATING TAG.**



DANGER



**DO NOT use this warming appliance
in the presence of flammable
anesthetic mixture (with air or
with oxygen or nitrous oxide).
THIS COULD CAUSE AN EXPLOSION!**

(Not category AP or APG equipment)

GENERAL WARNINGS



The unit requires special precautions regarding EMC (Electromagnetic Compatibility) and needs to be installed and put into service according to the EMC information provided in the accompanying documents.

Portable and mobile RF communications equipment can affect medical electrical equipment.

A risk of increased emissions or decreased immunity may result if the power cord attached is altered or a manufacturer supplied power cable is not used.

The unit should not be used adjacent to or stacked with other equipment.

Observe to verify normal operation if it is necessary to use adjacent to or stacked with other equipment.

Guidance and manufacturer's declaration – electromagnetic emissions

The units are intended for use in the electromagnetic environment specified below. The customer or the end user of this unit should assure that it is used in such an environment.

Emissions test	Compliance	Electromagnetic environment - guidance
RF emissions CISPR 11	Group 1	The unit uses RF energy only for internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class B	The unit is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/Flicker emissions IEC 61000-3-3	Complies	

Guidance and manufacturer's declaration – electromagnetic immunity

The unit is intended for use in the electromagnetic environment specified below. The customer or the end user of this unit should assure that it is used in such an environment.


Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Electromagnetic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	+2 kV for power supply lines	Main power quality should be that of a typical commercial or hospital environment. The unit does not have any input/output lines.
Surge IEC 61000-4-5	±1 kV differential mode ±2 kV common mode	±1 kV differential mode ±2 kV common mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5 % UT (>95 % dip in UT) for 0.5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles <5 % UT (>95 % dip in UT) for 5 sec	<5 % UT (>95 % dip in UT) for 0.5 cycle 40 % UT (60 % dip in UT) for 5 cycles 70 % UT (30 % dip in UT) for 25 cycles <5 % UT (>95 % dip in UT) for 5 sec	Mains power quality should be that of a typical commercial or hospital environment. If the user of the unit requires continued operation during power mains interruptions, it is recommended that the unit be powered from an uninterrupted power supply or a battery.
Power frequency (50/60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.

NOTE UT is the a.c. mains voltage prior to application of the test level.

The essential performance of the unit is to not exceed internal temperature of 150° F (+10%).

Guidance and manufacturer's declaration - electromagnetic emissions

The unit is intended for use in the electromagnetic environment specified below. The customer or the end user of this unit should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment - guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz	3 V/m	Portable and mobile RF communications equipment should be used no closer to any part of the unit, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d = [3.5/3] \sqrt{P}$ $d = [3.5/3] \sqrt{P}$ 80 MHz to 800 MHz $d = [7/3] \sqrt{P}$ 800 MHz to 2.5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	

NOTE 1 At 80 MHz and 800 MHz, the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the unit is used exceeds the applicable RF compliance level above, the unit should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the unit.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than [VI] V/m.

Guidance and manufacturer's declaration – electromagnetic immunity recommended separation distance between portable and mobile RF communications equipment and this unit.

The unit is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the unit can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the unit as recommended below, according to the maximum output power of the communications equipment.


Rated maximum output power of transmitter W	Separation distance according to frequency of transmitter m		
	150 kHz to 80 MHz	80 MHz to 800 MHz	800 MHz to 2.5 GHz
	$d = \left[\frac{3.5}{3} \right] \sqrt{P}$	$d = \left[\frac{3.5}{3} \right] \sqrt{P}$	$d = \left[\frac{7}{3} \right] \sqrt{P}$
0.01	0.117	0.117	0.233
0.1	0.369	0.369	0.738
1	1.167	1.167	2.333
10	3.689	3.689	7.379
100	11.667	11.667	23.333

For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.


NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

BLANKET CONTROL FEATURES



CONTROL PANEL BUTTONS



ON/OFF BUTTON

-  Press the ON/OFF button to power on the control.
- Press and hold the ON/OFF button for three (3) seconds to power off the control.


INTERIOR LIGHT BUTTON

-  Press INTERIOR LIGHT button to toggle blue interior LED light intensity to high, low, or off.

UP AND DOWN ARROW BUTTONS

-  Used to increase or decrease the temperature set-point. On units with optional timer installed the arrow buttons are used to increase or decrease time, date, auto-start, and auto-stop times.
- 

TEMPERATURE RECALL BUTTON



-  Press the TEMPERATURE RECALL button to view the actual cavity temperature. The display will show the actual cavity temperature for five (5) seconds before reverting back to displaying the current temperature set-point.

L.E.D. DIGITAL DISPLAY

The control has a four-digit L.E.D. display.


AUDIBLE BUTTON FUNCTION

The warmer's audible button function can be turned ON or OFF.


-  While the warmer is OFF, press and hold the DOWN ARROW button for four (4) seconds.
-  The display will show the current audible button status, 0 (OFF) or 1 (ON). Press the UP or DOWN arrows to toggle between the two states.

L.E.D. DISPLAY STATUS INDICATORS


OVERTEMP INDICATOR LIGHT

-  This indicator will illuminate and an alarm will sound if the control senses a temperature of 195°F (91°C) or higher. The alarm can be muted by pressing the ON/OFF button.

ERROR ACKNOWLEDGEMENT

-  To clear or acknowledge an error, press the ON/OFF button. Press the ON/OFF button to acknowledge the periodic alarm. If the alarm continues or returns, the warmer is still experiencing an error and may need service.

BUTTON LOCKOUT INDICATOR LIGHT


-  The lock indicator light will illuminate when the control lock feature is engaged. Press the ON/OFF button and UP ARROW simultaneously to lock the control. Press the ON/OFF button and DOWN ARROW button simultaneously to unlock the control.

POWER FAIL DETECTION


If the power fails for any reason while heating, the warmer will retain its current operating state in memory. The control will resume operating after the power is restored. In order to alert the user that the power has failed, the decimal place on the first digit will flash. Press the ON/OFF key to acknowledge the power failure and restoration.

NOTE: If the timer option is installed, the warmer must be off for more than 60 seconds to signal a power failure alarm. When acknowledging a power interruption, the display will show the length of time in hours and minutes that the control has been off due to the power outage.

TEMPERATURE FORMAT SELECTION

-  While the controller is in OFF state, press and hold the TEMPERATURE RECALL button key for four (4) seconds to display the current temperature scale. Press either the UP or DOWN arrows to toggle between Fahrenheit and Celsius.

OPTIONAL BLANKET CONTROL FEATURES



NOTE: If your warmer is not equipped with this optional automated timer feature or you choose not to use it, you do not need to set the time and this section can be skipped.

Unit time must be manually reset for daylight saving time.

Times are always displayed in 24-hour format (HH:MM). Midnight is 00:00. Noon is 12:00. 1:00 p.m. is 13:00.

If the display has changed due to user interaction, it will reset after five (5) seconds of inactivity.

If the start and stop times are the same, the warmer will recognize the OFF time only and the warmer will not turn on without user intervention. This is the best way to set the warmer for the days when it is not needed.

TIMER CONTROL PANEL

Blanket warmers with this option have four (4) timer buttons:

TIMER ON/OFF BUTTON



Press the TIMER ON/OFF button to initiate automatic start/stop operation mode. The ON/OFF indicator light next to this button will illuminate when mode is turned on.

UNIT TIME BUTTON



Press to view current time, date, and day and initiate changes to settings.

SET OR VIEW START TIME BUTTON



Press to view current automatic start time and initiate changes to set time.

SET OR VIEW STOP TIME BUTTON



Press to view current automatic stop time and initiate changes to set time.

RESET UNIT TIME BUTTON

Timer option must be OFF to reset the unit time.

NOTE: "E-60" is displayed if the clock is not set or the control has been off too long and the memory has corrupted.

1. If "E-60" is displayed on the screen, press either the UNIT TIME button or the ON/OFF button to acknowledge error.
2. Press and hold the UNIT TIME button for four (4) seconds until the auto-timer ON/OFF indicator light flashes slowly and the display shows the current set time. Press the UP or DOWN buttons to adjust the time by minute or press and hold the UP or DOWN button to adjust time more quickly.
3. Press the UNIT TIME button again to view or set the year. Press the UP or DOWN arrow keys to adjust the year. The year will always adjust by one (1).



4. Press the UNIT TIME button again to set the date. Press the UP or DOWN arrow keys to adjust the date.



5. Press the UNIT TIME button again to set the day code (d1-d7). Press the UP or DOWN arrows keys to adjust the day code. This is an optional step unless you are going to configure the auto-timer to start and/or stop at a different time for each day of the week. Typically Monday=d1, Tuesday=d2, etc., but this is user configurable.



SAME START AND STOP TIMES FOR THE WEEK

Timer option must be OFF to set start time.



1. Press the START TIME button to view the start time for that day.

2. Press and hold the START TIME button for four (4) seconds to set the current day's start time as the same time for every day of the week. The ON/OFF indicator light will flash slowly.



3. Press the UP or DOWN ARROW buttons to change the default start time minutes. Press and hold the arrows to adjust by hours.



4. Press the STOP TIME button to view the stop time for that day.

5. Press and hold the STOP TIME button for four (4) seconds to set the current day's stop time as the same time for every day of the week. The ON/OFF indicator light will flash slowly.



6. Press UP or DOWN ARROW buttons to change stop time. Press and hold the arrows to adjust by hours.



7. When both times are set, allow the ON/OFF indicator light to extinguish and then press and hold the Timer ON/OFF button until the ON/OFF indicator light stays on steadily.

Caution: Activating this mode overrides all individual day programming.

DIFFERENT START AND STOP TIMES FOR EACH DAY OF THE WEEK



Press and hold the UNIT TIME button for eight (8) seconds. The ON/OFF indicator light will flash rapidly. The display will show the current day. Adjust the day by pressing the UP or DOWN arrow keys.



1. Press the START TIME or STOP TIME button to display the start or stop times for that day. Press the UP or DOWN arrow keys to adjust the automatic start and stop times.



2. Press the UNIT TIME button again to display the day. Repeat steps 1-2 for each day of the week.



3. When all days and times are set, allow the ON/OFF indicator light to extinguish and then press and hold the Timer ON/OFF button until the ON/OFF indicator light stays on steadily.




4. The warmer will now turn on and off automatically at the times set for each day.

BLANKET CHAMBER OPERATION PROCEDURES

1. A DC Series appliance should be plugged into a hospital grade, 120V receptacle.

- DC150, DC250, DC350, DC400:
NEMA 5-15P; 15A - 125V Plug
- DC750: NEMA 5-20P; 20A - 125V Plug



2. **TURN ON THE POWER SWITCH**, which is located at the back of the appliance just above power cord.

3. **ACTIVATE CONTROL BY PRESSING THE ON/OFF KEY ON CONTROL PANEL ONCE.**

The digital L.E.D. display will indicate last temperature set-point of compartment.

4. **SET DESIRED TEMPERATURE.** To set the blanket warming temperature, press and hold the UP or DOWN ARROW keys to change the value shown in the display. The temperature set-point range is 37 - 82°C (98 - 180°F).

5. **LOAD THE CHAMBER WITH DRY, COTTON BLANKETS. DO NOT WARM ITEMS CONTAINING PLASTIC, RUBBER OR METAL SNAPS, STUDS, HOOKS, ETC.** Check that the epoxy-coated blanket insert is in place. This blanket insert **MUST** be used to hold blankets. Make certain the cabinet door is securely closed during use.

6. **ROTATE LOAD OF BLANKETS DAILY.**

Rotate the blankets at the bottom of the load to the top to ensure equal usage. Failure to rotate blankets can cause blankets to discolor.

Note: Avoid using flammable cabinet cleaning agents, as well as blanket cleaning agents that cause fabric to become brittle over time.

CAUTION

BLANKET SUPPORT INSERTS MUST BE USED WHEN WARMER IS IN USE.

CAUTION

**DO NOT OVERLOAD CABINET.
ALLOW 1" (25MM) GAP BETWEEN TOP INTERIOR OF UNIT OR SHELF AND BLANKETS.**



CLEANING AND PREVENTIVE MAINTENANCE

PROTECTING STAINLESS STEEL, EPOXY COATED AND PLASTIC SURFACES



It is important to guard against corrosion in the care of stainless steel surfaces. Harsh, corrosive, or inappropriate chemicals can completely destroy the protective surface layer of stainless steel, epoxy or plastic. Abrasive pads, steel wool, or metal implements will abrade surfaces causing damage to this protective coating and will eventually result in areas of corrosion. Even water, particularly hard water that contains high to moderate concentrations of chloride, will cause oxidation and pitting that result in rust and corrosion. In addition, many acidic spills left to remain on metal surfaces are contributing factors that will corrode surfaces.

Proper cleaning agents, materials, and methods are vital to maintaining the appearance and life of this appliance. Spilled items should be removed and the area wiped as soon as possible but at the very least, a minimum of once a day. Always thoroughly rinse surfaces after using a cleaning agent and wipe standing water as quickly as possible after rinsing.

ANNUAL PREVENTATIVE MAINTENANCE

1. Ensure that the correct Operation and Care Manual is available to all users.
2. Ensure that all users have been properly trained in unit's operation.
3. Do not overload cabinet.
 - **Blanket Warmer:** 1" (25mm) from top interior of unit
 - **Fluid Warmer:** See electrical/capacity page
4. Inspect condition of plug and cord. Replace if damaged.
5. Clean dust from outer vents surrounding the unit and around top of bonnet (if applicable).
6. Check door gasket integrity. Are there any tears? Is the gasket worn or loose? Make sure seal is tight to unit body. Replace gasket if integrity is compromised.
7. Check air temperature sensor mount on the interior of chamber. Is the guard in place? Are the wires in good condition?
8. Check insert assembly (depends on unit):
 - **Blanket Warmer:** Check the blanket support assembly and shelf. Is the assembly in place? Are any pieces missing?
 - **Fluid Warmer:** Check basket and side rail condition. Do baskets move smoothly and freely?
9. Check condition of casters or feet. Ensure components are secure and tightly threaded.
10. Check control panel overlay condition. Are there any tears or excessive wear on the graphic? Does the control work properly when buttons are pushed?
11. Check that all control and interior LEDs light up.
12. Is the set temperature comparable to the actual temperature displayed?

Contact service for immediate repair if any of the above problems exist.




CLEANING AGENTS

Use non-abrasive cleaning products designed for use on stainless steel surfaces. Cleaning agents must be chloride-free compounds and must not contain quaternary salts. Never use hydrochloric acid (muriatic acid) on stainless steel surfaces. Always use the proper cleaning agent at the manufacturer's recommended strength. Contact your local cleaning supplier for product recommendations.

CLEANING MATERIALS

The cleaning function can usually be accomplished with the proper cleaning agent and a soft, clean cloth. When more aggressive methods must be employed, use a non-abrasive scouring pad on difficult areas and make certain to scrub with the visible grain of surface metal to avoid surface scratches. Never use wire brushes, metal scouring pads, or scrapers to remove residue.



CAUTION	
  	TO PROTECT SURFACES, COMPLETELY AVOID THE USE OF ABRASIVE CLEANING COMPOUNDS, CHLORIDE BASED CLEANERS, OR CLEANERS CONTAINING QUATERNARY SALTS. NEVER USE HYDROCHLORIC ACID (MURIATIC ACID) ON STAINLESS STEEL. NEVER USE WIRE BRUSHES, METAL SCOURING PADS OR SCRAPERS.

CARE AND CLEANING

The cleanliness and appearance of this equipment will contribute considerably to its operating efficiency. Make certain the cabinet and door gasket are kept free of any debris that may accumulate. Good equipment that is kept clean works better and lasts longer.



CLEAN THE UNIT REGULARLY:



1. Disconnect the cabinet from the power source.
2. Remove all detachable items such as blanket insert. Clean these items separately.

NOTE: Avoid the use of abrasive cleaning compounds, chloride based cleaners, or cleaners containing quaternary salts. Never use hydrochloric acid (muriatic acid) on surfaces.



3. Clean the interior metal surfaces of the cabinet with a damp cloth and any mild commercial detergent. Avoid the use of abrasive cleaning compounds. Rinse surfaces by wiping with sponge and clean warm water. Remove excess water with sponge and wipe dry with a clean cloth or air dry. Leave doors open until interior is completely dry.
4. Interior can be wiped with a sanitizing solution after cleaning and rinsing. This solution must be approved for use on stainless steel surfaces. Replace blanket support assembly.
5. Clean the exterior of the cabinet with a cleaner recommended for powder coated surfaces.
6. Clean the window glass with a standard commercial glass cleaner.
7. Wipe control panel, door vents, door buttons, and door gaskets thoroughly since these areas can harbor debris.
8. Wipe door gaskets and control panel dry with a clean, soft cloth.

Always follow appropriate state or local health (hygiene) regulations regarding all applicable cleaning and sanitation requirements.



(IPX-0 - Listed as Ordinary)

TROUBLESHOOTING GUIDE

NOTE: If your unit is not operating properly, check the following before calling your authorized service agent. Check the power applied to the unit. Verify female end of plug is securely seated in unit and that the male end of plug is in an appropriate, functioning outlet. Check fuses. (See "Fuse Replacement" section in manual.) Has the high limit manual reset tripped? If so, reset. (See "Manual Reset Instructions" below.)

Do not attempt to repair or service beyond this point. Contact manufacturer for nearest authorized service agent. Repairs made by any other service agent without prior authorization by manufacturer will void the warranty on the unit.

This chart is provided for the assistance of qualified technicians only and is not intended for use by untrained or unauthorized service personnel.

Code	Refers to	Action Required
E-10 ES10 ES20 ES30 ES40 ES50 ES60 ES70	Cavity sensor Pad sensor 1 Pad sensor 2 Pad sensor 3 Pad sensor 4 Pad sensor 5 Pad sensor 6 Pad sensor 7	Sensor is shorted. Software disengages heating pads. User must acknowledge error by pressing ON/OFF button. If error persists, a qualified service technician should test sensor. <ul style="list-style-type: none"> To test sensor: Detach the sensor from unit. Use an Ohm meter to measure the resistance of the sensor. Check sensor at 0°C (32°F) using a container of ice water. If reading is 32.6 KOhm ±1.5KOhm, replace display. If reading is ±2KOhm, replace sensor. Check wires for integrity. Check for proper and secure connections at the control and terminal block. If necessary, re-secure the faulty connections. Call Service if error persists.
E-11 ES11 ES21 ES31 ES41 ES51 ES61 ES71	Cavity sensor Pad sensor 1 Pad sensor 2 Pad sensor 3 Pad sensor 4 Pad sensor 5 Pad sensor 6 Pad sensor 7	Sensor is open. Software disengages heating pads. User must acknowledge error by pressing ON/OFF button. If error persists, a qualified service technician should test sensor. <ul style="list-style-type: none"> To test sensor: Detach the sensor from unit. Use an Ohm meter to measure the resistance of the sensor. Check sensor at 0°C (32°F) using a container of ice water. If reading is 32.6 KOhm ±1.5KOhm, replace display. If reading is ±2KOhm, replace sensor. Check wires for integrity. Check for proper and secure connections at the control and terminal block. If necessary, re-secure the faulty connections. Call Service if error persists.
P130 P230 P330 P430 P530 P630 P730	Pad 1 Pad 2 Pad 3 Pad 4 Pad 5 Pad 6 Pad 7	Heating pad has not reached set-point temperature within one hour. User must acknowledge error by pressing ON/OFF button. If error persists, a qualified service technician should test the pad(s). <ul style="list-style-type: none"> Turn unit OFF and unplug it from AC power. Use Ohm meter to measure resistance between L:(Line) & N:(Neutral) leads of heater pad. The Ohm readings shall be: <ul style="list-style-type: none"> 120V = 72 Ohm ±10% 230V = 288 Ohm ±10% Call Service if error persists.
E-31	Cavity sensor	Sensor reading is above maximum allowable temperature set-point. Call Service.
P131 P231 P331 P431 P531 P631 P731	Pad sensor 1 Pad sensor 2 Pad sensor 3 Pad sensor 4 Pad sensor 5 Pad sensor 6 Pad sensor 7	<ul style="list-style-type: none"> Software disengages heating pads. User must acknowledge error by pressing ON/OFF button. Allow unit to cool. Call Service if error persists.
E-60	Real Time Clock Checksum Error	Real Time Clock rechargeable battery backup has discharged. <ul style="list-style-type: none"> Plug unit into outlet for 30 minutes. See "Timer Control Panel" section in manual to reset clock.
*E-61	Real Time Clock	Real Time Clock not responding. Call Service if error persists.
E-62	Real Time Clock	Timer overlay is present, but no real time clock is detected. Call Service.
*E-70	Pad Count Error	More pads detected than set for. Hold ON/OFF button for 12 seconds until display shows "Pad#" (# = number of pads selected [3-7]). Press UP or DOWN arrow to adjust to correct number of pads. (DC150 & DC250 = 3 pads, DC350 & DC400 = 4 pads, DC750 = 7 pads)
*E-71	Personality Error/ Incorrect Dip Switch Setting	Call Service
E-80	EEPROM Error	Call Service.
E-83	EEPROM Error	Call Service for help resetting the control.
E-87	EEPROM Error	Stored offsets corrupted. Offsets reset to 0. Control may need a recalibration. Possible bad EEPROM. Call Service if error persists.
E-90	Button stuck	A button has been held down for >60 seconds. Adjust control. Error will reset when the problem has been resolved.

Note: All non-critical codes can be cleared using the ON/OFF button. Critical errors (marked with a *) can only be cleared by turning the power switch on the back of the unit off.

Manual Reset Instructions: After allowing unit to cool, remove three (3) screws from top cover (on upper edge of the back of unit). Slide top panel toward back of unit (approximately 1/2" [13mm]) and then lift up. Locate the two (2) manual reset buttons in middle of the bonnet. (See "MANUAL RESETS" on the service views on following pages.) Firmly push reset button(s). You will hear an audible click when the buttons are reset. If reset button trips again while unit is running, contact a qualified service technician.

SERVICE

FULL ASSEMBLY


LOC	DESCRIPTION	DC150 P/N	QTY	DC250 P/N	QTY	DC350 P/N	QTY	DC400 P/N	QTY	DC750 P/N	QTY
1.	STAND-OFFS/SPACERS	SP-33849	11	SP-33849	11	SP-33849	11	SP-33849	11	SP-33849	11
2.	CONTROL BOARD	STANDARD BA-34867	1	BA-34867	1	BA-34867	1	BA-34867	1	BA-34867	1
*	W/OPTIONAL TIMER CLOCK	BA-34869	1	BA-34869	1	BA-34869	1	BA-34869	1	BA-34869	1
3.	RELAY BOARD	3-ZONE BA-34950	1	BA-34950	1	—	—	—	—	—	—
*	4-ZONE	—	—	—	—	BA-34951	1	BA-34951	1	—	—
*	7-ZONE	—	—	—	—	—	—	—	—	BA-34954	1
4.	POWER SUPPLY BOARD	BA-34965	1	BA-34965	1	BA-34965	1	BA-34965	1	BA-34965	1
5.	TOP COVER	1011029	1	1011266	1	1011266	1	1011475	1	1011475	1
6.	POWER SWITCH ASSEMBLY	SW-34911	1	SW-34911	1	SW-34911	1	SW-34911	1	SW-34769	1
7.	FUSES (F1, F2) FAST ACTING	250V, 10A FU-34955	2	FU-34955	2	FU-34955	2	FU-34955	2	—	—
	250V, 15A	—	—	—	—	—	—	—	—	FU-35094	2
8.	FUSE HOLDER	—	—	—	—	—	—	—	—	FU-35033	2
9.	INLET, POWER	—	—	—	—	—	—	—	—	IT-29599	1
10.	PLUG, 1/2" ROUND	PG-3344	1	PG-3344	1	PG-3344	1	PG-3344	1	PG-3344	1
11.	POWER CORD	CD-35030	1	CD-35030	1	CD-35030	1	CD-35030	1	CD-34993	1
12.	REAR COVER	5013520	1	5013521	1	5013522	1	5013643	1	5013801	1
13.	SCREWS, HINGES (UPPER AND LOWER)	SC-27046	6	SC-27046	6	SC-27046	6	SC-27046	6	SC-27046	6
14.	HINGE ASSEMBLY, UPPER	5015522	1	5015522	1	5015522	1	5015522	1	5015522	1
15.	SPRING, DOOR	SD-29340	1	SD-29467	1	SD-29467	1	SD-29467	1	SD-29467	1
16.	SIDE COVER	1011024	2	1011269	2	1011288	2	1011472	2	1011588	2
17.	CLIP, HEATING PAD	CL-22259	18	CL-22259	18	CL-22259	18	CL-22259	24	CL-22259	36
18.	INSULATION ASSEMBLY	5014403	1	5014404	1	5014404	1	5015511	1	5015511	1
19.	HEATING PAD ASSEMBLY INCLUDES HEAT PAD ELEMENT (EL-34838) AND REQUIRED CONNECTORS	5013210	3	5013210	3	5013210	4	5013210	4	5013210	7
20.	HEAT PAD SENSOR	PR-34837	3	PR-34837	3	PR-34837	4	PR-34837	4	PR-34837	7
21.	THERMOSTAT, AUTO RESET	TT-34908	3	TT-34908	3	TT-34908	4	TT-34908	4	TT-34908	7
22.	BUMPER FEET, RUBBER	FE-29203	4	FE-29203	4	FE-29203	4	FE-29203	4	FE-29203	4
*	DOOR ASSEMBLY	5011206	1	5011584	1	5011606	1	5012014	1	5012248	1
23.	DOOR GASKET	GS-29188	1	GS-29189	1	GS-29190	1	GS-29191	1	GS-29344	1
24.	HINGE ASSEMBLY, LOWER	5015523	1	5015523	1	5015523	1	5015523	1	5015523	1
25.	BOTTOM COVER	5013523	1	5013524	1	5013524	1	5013644	1	5013644	1
26.	BLANKET SUPPORT, BOTTOM	1011045	1	1011265	1	1011265	1	1011471	1	1011471	1
27.	BLANKET SUPPORT, SIDES	1011046	2	1012866	2	1011358	4	1011470	2	1011587	4
28.	HEATED SHELF	—	—	—	—	5015540	1	—	—	5015545	1
29.	SHELF, PANEL COVER	—	—	—	—	1011230	1	—	—	1011465	1
30.	LED ASSEMBLY	5012953	1	5012953	1	5012953	2	5012953	1	5012953	2
31.	CLIP, PANEL RETAINER	CL-29193	6	CL-29193	6	CL-29193	6	CL-29193	6	CL-29193	6
32.	INTERFACE ASSEMBLY†										
	18" (457mm); STANDARD	5012323	1	5012323	1	5012323	1	—	—	—	—
*	18" (457mm) w/TIMER; OPTIONAL	5012690	1	5012690	1	5012690	1	—	—	—	—
*	24" (610mm); STANDARD	—	—	—	—	—	—	5012325	1	5012325	1
*	24" (610mm) w/TIMER; OPTIONAL	—	—	—	—	—	—	5012691	1	5012691	1
33.	CLIPS, CHASSIS	CL-35031	4	CL-35031	4	CL-35031	4	CL-35031	4	CL-35031	4
34.	THERMOSTAT, MANUAL RESET	TT-34907	2	TT-34907	2	TT-34907	2	TT-34907	2	TT-34907	2
35.*	SCREWS, COVERS & SHELF	SC-22271	**	SC-22271	**	SC-22271	**	SC-22271	**	SC-22271	**
36.*	WIRE DIAGRAM	77321	1	77321	1	77321	1	77321	1	77321	1

*NOT SHOWN ** VARIES † SEE PG 16 FOR EXPLODED VIEW OF INTERFACE

PART NUMBERS AND DRAWINGS ARE SUBJECT TO CHANGE WITHOUT NOTICE.

SERVICE

FULL ASSEMBLY (DC150 SHOWN)




PART NUMBERS AND DRAWINGS ARE SUBJECT TO CHANGE WITHOUT NOTICE.

SERVICE

INTERFACE ASSEMBLY (18" WITH TIMER OPTION SHOWN)

LOC	DESCRIPTION	P/N	QTY	
1.	PANEL, INTERFACE MOUNT	18" (DC150, DC250 & DC350)	1011336	1
		24" (DC400 & DC750)	1011482	1
2.	SPACER, SPEAKER BOARD		SP-29425	3
3.	SPEAKER BOARD		BA-34979	1
4.	SCREW, SPEAKER BOARD, M3x0.5x6mm		SC-22270	3
5.	SCREWS, M4x0.7x12mm		SC-29333	15
6.	WASHER		WS-22323	1
7.	WASHER		WS-25056	1
8.	PANEL, FRONT TOP TRIM	18" (DC150, DC250 & DC350)	1011032	1
		24" (DC400 & DC750)	1011481	1
9.	WASHER		WS-23991	1
10.	SPRING, COMPRESSION		SD-29371	1
11.	CONTROL ASSEMBLY SERVICE KIT (INCLUDES COVERS, OVERLAY & SCREWS)		5015524	1
12.	CIRCUIT BOARD		BA-34868	1
13.	BUTTON RETAINER PLATE		1011149	1
14.	COVER, BUTTON W/ TIMER (OPTIONAL)	18" (DC150, DC250 & DC350)	CV-29434	1
*		24" (DC400 & DC750)	CV-35433	1
*	COVER, BUTTON	18" (DC150, DC250 & DC350)	CV-29256	1
*		24" (DC400 & DC750)	CV-29345	1
15.	LATCH, BUTTON		LT-29176	1
16.	OVERLAY, TIMER (OPTIONAL)		PE-34864	1

*NOT SHOWN




PART NUMBERS AND DRAWINGS ARE SUBJECT TO CHANGE WITHOUT NOTICE.

SERVICE

LATCH MECHANISMS & CAVITY SENSOR


LOC	DESCRIPTION	DC150 P/N	QTY	DC250 P/N	QTY	DC350 P/N	QTY	DC400 P/N	QTY	DC750 P/N	QTY
1.	DOOR BUTTON LATCH PIVOT ASSEMBLY	5013640	1	5013640	1	5013640	1	5013640	1	5013640	1
2.	LED ASSEMBLY	5012953	1	5012953	1	5012953	2	5012953	1	5012953	2
3.	BUSHING	BU-3006	1	BU-3006	1	BU-3006	1	BU-3006	1	BU-3006	1
4.	SENSOR, CAVITY	SN-34976	1	SN-34976	1	SN-34976	1	SN-34976	1	SN-34976	1
5.	WASHER, FLAT	WS-2420	3	WS-2420	3	WS-2420	3	WS-2420	3	WS-2420	3
6.	BRACKET, SENSOR MOUNT	1012011	1	1012011	1	1012011	1	1012011	1	1012011	1
7.	SCREWS, SENSOR BRACKET, M4x0.7x6mm	SC-22271	2	SC-22271	2	SC-22271	2	SC-22271	2	SC-22271	2
8.	GUIDE, LOW PROFILE RAIL	GI-29204	1	GI-29204	1	GI-29204	1	GI-29204	1	GI-29204	1
9.	GUIDE, CARRIAGE	GI-29205	2	GI-29205	2	GI-29205	2	GI-29205	2	GI-29205	2
10.	SCREWS, M4x0.7x6mm	SC-22273	8	SC-22273	8	SC-22273	8	SC-22273	8	SC-22273	8
11.	LATCH	LT-29174	2	LT-29174	2	LT-29174	2	LT-29174	2	LT-29174	2
12.	PLATE, LATCH	1011197	1	1011263	1	1011295	1	1011263	1	1011850	1
13.	HOUSING, LATCH SPRING SUPPORT	1011163	1	1011163	1	1011163	1	1011163	1	1011163	1
14.	SPRING	SD-28513	1	SD-28513	1	SD-28513	1	SD-28513	1	SD-28513	1
15.	BUSHING	BU-29206	3	BU-29206	3	BU-29206	3	BU-29206	3	BU-29206	3
16.	CLIP	CL-29257	3	CL-29257	3	CL-29257	3	CL-29257	3	CL-29257	3




PART NUMBERS AND DRAWINGS ARE SUBJECT TO CHANGE WITHOUT NOTICE.

ZONE HEATING PAD LOCATIONS

DC150, DC250, DC350 Heating Pad Locations




DC400, DC750 Heating Pad Locations



FUSE REPLACEMENT

Fuse replacement for DC150, DC250, DC350 & DC400 models:

1. Unplug power cord from wall outlet and from power switch assembly.
2. Fuse compartment is located directly below the plug receptacle. Use fingernail or thin implement to flip compartment door open.
3. Use fingernail or thin implement to pull fuse drawer out from compartment.
4. Use a thin implement to push fuses up out of drawer.
5. Replace with new fuse.
6. Push drawer back into compartment.
7. Close compartment door.



⚠ WARNING

FOR PROTECTION AGAINST FIRE AND ELECTRICAL SHOCK.


USE ONLY UL LISTED 10A, 250V FAST ACTING FUSES, 5MM X 20MM (F1, F2).

USE ONLY HOSPITAL GRADE CORD. SEE OPERATION AND CARE MANUAL FOR LIST OF AVAILABLE CORDS.

LA-29452

Fuse replacement for DC750 models:

1. Unplug power cord from wall outlet and power inlet on unit.
2. Fuse compartments are located next to the plug receptacle. Unscrew holder cap.
3. Pull fuses out of holder.
4. Replace with new fuses.
5. Screw cap back on.



⚠ WARNING

FOR PROTECTION AGAINST FIRE AND ELECTRICAL SHOCK.

USE ONLY UL LISTED 15A, 250V FAST ACTING FUSES, 6.3MM X 32MM (F1, F2).

USE ONLY HOSPITAL GRADE CORD. SEE OPERATION AND CARE MANUAL FOR LIST OF AVAILABLE CORDS.

LA-35440


OPTIONS & ACCESSORIES PARTS LIST*

*NOT SHOWN

OPTIONS & ACCESSORIES	P/N	QTY
STEM CASTER, SWIVEL 3" (76mm)	CS-20500	2
STEM CASTER, BRAKE 3" (76mm)	CS-29499	2

PART NUMBERS AND DRAWINGS ARE SUBJECT TO CHANGE WITHOUT NOTICE.

Refer to wire diagram under top cover of unit for most current version



REV	ECO	DESCRIPTION	DATE	APP
13	EC0265	UPDATE SWITCH NUMBERS	1/3/2013	JJD
12	EC0844	UPDATE WIRE NUMBERS	8/20/2012	JJD
11	EC0224	UPDATE RB PART NUMBERS	07/16/2012	JJD

DELUXE	WARMER_120/220/230V_60/50Hz
SHEET	DATE
1_OF_1	DATE

BY: PKR DATE: 01/27/10 DWG: 77321

NOTE #1: SEE 5012058, 5014425, AND 5013351 FOR COMMON WIRE ASSEMBLIES

CAPACITY	WATTS							AMPS							
	1.5F ³	2.5F ³	3.5F ³	4.0F ³	7.5F ³	1.5F ³	2.5F ³	3.5F ³	4.0F ³	7.5F ³	1.5F ³	2.5F ³	3.5F ³	4.0F ³	7.5F ³
DESIGN	120V	615W	815W	815W	1415W	5.1A	5.1A	5.1A	6.8A	6.8A	11.8A	11.8A	11.8A	11.8A	11.8A
	220V	540W	540W	715W	1240W	2.5A	2.5A	2.5A	3.3A	3.3A	5.6A	5.6A	5.6A	5.6A	5.6A
	230V	615W	615W	815W	1415W	2.7A	2.7A	2.7A	3.6A	3.6A	6.2A	6.2A	6.2A	6.2A	6.2A
	120V	600W	600W	800W	1400W	5.0A	5.0A	5.0A	6.7A	6.7A	11.7A	11.7A	11.7A	11.7A	11.7A
	220V	550W	550W	700W	1250W	2.5A	2.5A	2.5A	3.2A	3.2A	5.7A	5.7A	5.7A	5.7A	5.7A
	230V	600W	600W	800W	1400W	2.6A	2.6A	2.6A	3.5A	3.5A	6.1A	6.1A	6.1A	6.1A	6.1A

TRANSPORTATION DAMAGE AND CLAIMS



All Enthermics Medical Systems equipment is sold F.O.B. shipping point, and when accepted by the carrier, such shipments become the property of the consignee.

Should damage occur in shipment, it is a matter between the carrier and the consignee. In such cases, the carrier is assumed to be responsible for the safe delivery of the merchandise, unless negligence can be established on the part of the shipper.

1. Make an immediate inspection while the equipment is still in the truck or immediately after it is moved to the receiving area. Do not wait until after the material is moved to a storage area.
2. Do not sign a delivery receipt or a freight bill until you have made a proper count and inspection of all merchandise received.
3. Note all damage to packages directly on the carrier's delivery receipt.
4. Make certain the driver signs this receipt. If he refuses to sign, make a notation of this refusal on the receipt.
5. If the driver refuses to allow inspection, write the following on the delivery receipt:

Driver refuses to allow inspection of containers for visible damage.

6. Telephone the carrier's office immediately upon finding damage, and request an inspection. Mail a written confirmation of the time, date, and the person called.
7. Save any packages and packing material for further inspection by the carrier.
8. Promptly file a written claim with the carrier and attach copies of all supporting paperwork.

We will continue our policy of assisting our customers in collecting claims which have been properly filed and actively pursued. We cannot, however, file any damage claims for you, assume the responsibility of any claims, or accept deductions in payment for such claims.

ENTHERMICS MEDICAL SYSTEMS LIMITED WARRANTY

Enthermics Medical Systems warrants to the original purchaser that any original part that is found to be defective in material or workmanship will, at our option, subject to provisions hereinafter stated, be replaced with a new or rebuilt part.

The labor warranty remains in effect one (1) year from installation or fifteen (15) months from the shipping date, whichever occurs first.

The original parts warranty for the cavity fan motor remains in effect one (1) year from installation of appliance or fifteen (15) months from the shipping date, whichever occurs first. The original parts warranty on all other parts remains in effect three (3) years from installation of appliance or thirty-nine (39) months from the shipping date, whichever occurs first.

This warranty does not apply to:

1. Calibration
2. Equipment damage caused by accident, shipping, improper installation or alteration.
3. Equipment used under conditions of abuse, misuse, carelessness or abnormal conditions including equipment subjected to harsh or inappropriate chemicals including but not limited to compounds containing chloride or quaternary salts, poor water quality, or equipment with missing or altered serial numbers.
4. Any losses or damage resulting from malfunction, including loss of contents or consequential or incidental damages of any kind.
5. Equipment modified in any manner from original model, substitution of parts other than factory authorized parts, removal of any parts including legs, or addition of any parts.
6. Collateral or incidental damage as a direct result of servicing equipment built into a wall structure is not covered under warranty. It is the responsibility of the owner to bear all expense related to structural repairs including, but not limited to, external electrical connections and wiring, and the removal or replacement of caulk, grout, tile, or wall covering of any kind. A service access panel for built-in equipment installations is strongly recommended.

This warranty is exclusive and is in lieu of all other warranties, expressed or implied, including the implied warranties of merchantability and fitness for purpose. In no event shall the Company be liable for loss of use, loss of revenue, or loss of contents or revenue, or for indirect or consequential damages. This warranty is in lieu of all other warranties expressed or implied and Enthermics Medical Systems neither assumes or authorizes any persons to assume for it any other obligation or liability in connection with Enthermics Medical Systems equipment.

Enthermics Medical Systems

Record the model and serial numbers of the unit for easy reference. Always refer to both model and serial numbers in your correspondence regarding the unit.

Model: _____
Serial Number: _____
Purchased From: _____
Date Installed: _____ Voltage: _____

Warranty Effective November 1, 2012



ENTHERMICS Medical Systems

Corporate Headquarters:

W164 N9221 Water Street • P.O. Box 443
Menomonee Falls, Wisconsin 53052-0443 • U.S.A.

PHONE: 262.251.8356
www.enthermics.com

800.T.O.B.WARM U.S.A./CANADA
800.862.9276 U.S.A./CANADA

FAX: 262.251.7067
800.329.8744